



Pre-MRI Session

Step 1: Register new patient on MRI Operator Computer

- New patient must be registered prior to launching **FIRMM**.

Step 2: FIRMM session start on MRI Operator Computer

1. After participant is registered, access the start menu (*Ctrl+Esc*).
2. Click FIRMM session start.
3. If a text file comes up, close it. Instead, right-click FIRMM session start and click **Open**.

Note 1: FIRMM session start only applies to current registered patient.

Note 2: Scanner will only send data from “applied” scans AFTER the FIRMM session start.

Step 3: Login to **FIRMM on FIRMM Computer**

1. Login as “FIRMM user”: username _____ password _____
2. Open a Terminal and type **FIRMM** then press Enter to bring up the **FIRMM** window.

Step 4: Start **FIRMM on FIRMM Computer**

1. Press **Start** on right.
2. **FIRMM** will ask for your current scan. If you do not see it, click **Refresh** to check again.

Note: You can only see the current scan after FIRMM session start on the scanner and after the first DICOMs have transferred. If you have not started scanning, you will not see it yet.

3. Select the session you want to monitor and click **Run**.

Post-MRI Session

Close **FIRMM**

- When the scanning session is over, close the **FIRMM** window on the FIRMM Computer.

Note: To run another MRI session with **FIRMM**, start over again from **Step 1** above.

FIRMM Motion Data: CSVs with FD and motion numbers are stored at ‘/home/firmmproc/FIRMM/outgoing/FIRMM_logs’ on the FIRMM computer.